Highlighting innovative design features and useful applications information for

Thermo Scientific Centrifuges and Rotors


smart notes

design and innovation









How does Auto-Lock rotor exchange facilitate rotor placement and improve safety, compared with traditional rotor tie-down systems?

In only 3 seconds, with just the push of a button, Thermo Scientific Auto-Lock rotor exchange allows you to install and remove rotors on select¹ Thermo Scientific superspeed or general purpose centrifuges. With Auto-Lock® rotor exchange, you can have the confidence that your rotor is safely and securely locked in place, and have the flexibility to change rotors and applications in the blink of an eye.

Unlike traditional rotor tie-down systems, with push-button Auto-Lock rotor exchange, rotors no longer need to be bolted onto a centrifuge motor shaft. As a result:

- 1 No tools are required
- 2] No complicated rotor changing procedure is necessary
- 3 Swap between protocols instantly for application flexibility
- 4 Clear chamber access in only seconds for quick and efficient cleaning

Auto-Lock is identified on Therma Scientific superspeed

Auto-Lock is identified on Thermo Scientific superspeed centrifuges with a **black push button** and on Thermo Scientific general purpose centrifuges with a **green push button**.

¹ Auto-Lock is available with the following Thermo Scientific centrifuge series:

Superspeed	General Purpose	
 Sorvall® LYNX 6000/4000 	 Sorvall Legend® X1/XT/XF Sorvall ST 16/40 	 Heraeus® Multifuge® X1/X3 Heraeus Megafuge® 16/40
		• SL 16/40

Rotors featuring Auto-Lock are not interchangeable between superspeed and general purpose benchtop platforms.

In today's busy laboratories, safe and efficient sample processing is essential to getting answers faster.

Traditional Rotor Placement











1 Place rotor


2 Manually tighten

3 Secure with tool

4] Confirm locked down

Auto-Lock Rotor Exchange





1 Place rotor with confidence that it is securely locked and will not loosen during a run

Push-button Security and Application Flexibility

In a shared laboratory setting with multiple users and a variety of processing requirements, or where bench space is at a premium, Auto-Lock rotor exchange provides the flexibility to quickly swap between your sample preparation protocols, from 0.2 mL through 1 Liter, all with the same versatile centrifuge. The easy and reliable push-button Auto-Lock rotor exchange system replaces the need for tools to bolt down the rotor, shortening run set-up time. Also, as the rotor automatically and securely locks itself to the centrifuge, the need for hand-tightening is eliminated, improving safety and confidence that the rotor will not loosen during a run.

Routine Cleaning Made Easy

Centrifuge chambers can harbor microscopic contamination from the environment and from spilled samples. To ensure a healthy working environment and prevent contamination, laboratories require that the centrifuge chamber is cleaned and disinfected regularly. Auto-Lock rotor exchange allows removal of the rotor in only 3 seconds, providing quick and easy access to the centrifuge chamber for the routine cleaning. And once cleaned, the rotor can be securely locked in place again in only 3 seconds.

Summary

Thermo Scientific Auto-Lock rotor exchange delivers secure, trouble-free rotor installation and removal, easy routine cleaning for safe sample processing and 3-second rotor exchange to simplify performance and evolve with the changing needs of your laboratory.

Experience Auto-Lock rotor exchange now at thermoscientific.com/auto-lock

© 2012 Thermo Fisher Scientific Inc. All rights reserved. HERAEUS is a registered trademark of Heraeus Holding GmbH licensed to Thermo Fisher Scientific. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Australia +61 39757 4300 **Austria** +43 1 801 40 0 **Belgium** +32 53 73 42 41

China +86 21 6865 4588 or +86 10 8419 **France** +33 2 2803 2180

Germany national toll free $0800\ 1\ 536\ 376$ Germany international $+49\ 6184\ 90\ 6000$

India toll free 1800 22 8374 India +91 22 6716 2200 Italy +32 02 95 05 92 54 Japan +81 45 453 9220 Netherlands +31 76 579 55 55 New Zealand +64 9 980 6700 Nordic/Baltic/CIS countries +358 9 329 10200 Russia +7 812 703 42 15 Spain/Portugal +34 93 223 09 18 Switzerland +41 44 454 12 22 UK/Ireland +44 870 609 9203 USA/Canada +1 866 984 3766

Other Asian countries +852 2885 4613 **Countries not listed** +49 6184 90 6000

thermoscientific.com/rotor



Part of Thermo Fisher Scientific